

IST Sample Preparation Kit (8 reactions) PROTOCOL - Urine

INTRODUCTION

Sample preparation is one of the essential steps of bottom-up proteomics. The PreOmics sample preparation kit is designed to assist you achieving best results with few sample preparation steps and little hands-on time. For sample-specific protocols and optimization contact us or visit our website at www.preomics.com.

KIT CONTENTS

The kit contains all you need to perform a complete sample preparation. It includes all chemicals to denature, reduce and alkylate proteins, as well as the enzymes to perform a tryptic digestion and a final cleanup.


| Component | Cap | Quantity | Buffer Properties | | | | Description | Storage |
|------------|----------|----------|-------------------|--------|-------|----------|--|---------|
| | | | Organic | Acidic | Basic | Volatile | | |
| DIGEST | red | 2x | | | | | Trypsin/LysC mix to digest proteins. | -20°C |
| RESUSPEND | yellow | 1x 2 mL | | | | ☉ | Protease reconstitution buffer for enzymes. | RT |
| LYSE | brown | 1x 1 mL | | | | ☉ | Denatures, reduces and alkylates proteins. | RT |
| STOP | black | 1x 1 mL | ☉ | ☉ | | ☉ | Stops the enzymatic activity. | RT |
| WASH 0 | no color | 1x 2 mL | ☉ | ☉ | | ☉ | Clean up peptides from tetrapyrrole contaminants. | RT |
| WASH 1 | blue | 1x 2 mL | ☉ | ☉ | | ☉ | Clean up peptides from hydrophobic contaminants. | RT |
| WASH 2 | green | 1x 2 mL | | ☉ | | ☉ | Clean up peptides from hydrophilic contaminants. | RT |
| ELUTE | violet | 1x 2 mL | ☉ | | ☉ | ☉ | Elutes the peptides from the cartridge. | RT |
| LC-LOAD | white | 1x 1 mL | | ☉ | | ☉ | Load peptides on reversed-phase LC-MS column. | RT |
| CARTRIDGES | | 8x | | | | | Cartridge for 1 to 100 µg protein starting material. | RT |
| WASTE | | 8x | | | | | Tube for collecting waste after washing steps. | RT |
| COLLECTION | | 8x | | | | | Tube for collecting peptides after elution. | RT |
| ADAPTER | | 8x | | | | | Enables a cartridge to be placed into a tube. | RT |
| CAP | | 8x | | | | | Cap to (optionally) close the cartridge's bottom. | RT |

PRE-REQUISITES

Common lab equipment is required for the sample preparation.

| EQUIPMENT | QUANTITY AND DESCRIPTION |
|--------------------|---|
| PIPETTE | Careful sample handling and pipetting reduces contaminations and improves quantification. |
| SAMPLE | 100 mL urine. |
| CENTRIFUGAL FILTER | 3kDa protein concentration ultrafiltration system (e.g. Merck Millipore ACS500302). |
| HEATING BLOCK | Two heating blocks are recommended to support protein denaturation and digestion. |
| CENTRIFUGE | Eppendorf tube centrifuges are necessary for loading, washing and elution. |
| SONICATOR | If you have DNA in your sample, shear it by sonication (e.g. Diagenode Bioruptor®). |
| SPEED-VAC | Vacuum manifolds evaporate volatile buffers from the eluate before LC-MS. |
| ULTRASONIC BATH | Optional: can be used to resuspend peptides. |

PROCEDURE



IST Sample Preparation Kit (8 reactions)

PROTOCOL - Urine



PROTOCOL

1. PREPARE AND LYSE

- 1.1. Load 10 mL urine onto a centrifugal filter (not provided) and concentrate to 1 mL. Repeat until 100 mL are loaded.
- 1.2. Concentrate the sample to 100 µL and add 900 µL LYSE ● to your sample. Mix gently.
- 1.3. Concentrate the solution to 50 µL, mix thoroughly and transfer the sample to a fresh tube.
- 1.4. Place the sample in a HEATING BLOCK (95°C; 1,000 rpm; 10 min). *NOTE 1*
- 1.5. Optional: Spin down droplets (RT; max. 300 rcf; 10 sec).
- 1.6. If you expect DNA, shear sample in a SONICATOR (10 cycles; 30 sec ON/OFF).
- 1.7. Use ADAPTER to place CARTRIDGE in WASTE tube. Label all tubes.
- 1.8. Transfer sample to CARTRIDGE and cool down (RT). Be careful not to damage the bottom layer of CARTRIDGE.

2. DIGEST

- 2.1. Add 210 µL RESUSPEND ● to DIGEST ● (1 tube for 4 reactions), shake (RT; 500 rpm; 10 min), pipette up/down.
- 2.2. Add 50 µL DIGEST ● to CARTRIDGE and place it in a pre-heated HEATING BLOCK (37 °C; 500 rpm; 1-3 hours).
- 2.3. Add 100 µL STOP ● to CARTRIDGE (precipitation may occur), shake (RT; 500 rpm; 1 min), pipette up/down. *SP*

3. PURIFY

- 3.1. Spin CARTRIDGE in a CENTRIFUGE (3,800 rcf; 1-3 min). If needed, adjust values to ensure complete flow-through.
- 3.2. Load 200 µL WASH 0 (no color), repeat step 3.1.
- 3.3. Add 200 µL WASH 1 ● to CARTRIDGE, repeat step 3.1.
- 3.4. Add 200 µL WASH 2 ● to CARTRIDGE, repeat step 3.1. *SP*
- 3.5. Use ADAPTER to place CARTRIDGE in a fresh COLLECTION tube. Label all tubes.
- 3.6. Add 100 µL ELUTE ● to CARTRIDGE, repeat step 3.1., keep flow-through in COLLECTION tube.
- 3.7. Repeat step 3.6., keep flow-through in the same COLLECTION tube.
- 3.8. Discard CARTRIDGE and place COLLECTION tube in a SPEED-VAC (45 °C; until completely dry).
- 3.9. Add LC-LOAD ○ to COLLECTION tube. Aim for 1 g/L concentration (e.g. 100 µL to 100 µg protein starting material).
- 3.10. Sonicate COLLECTION tube in an ULTRASONIC BATH (5 min) or shake (RT; 500 rpm; 5 min). *SP*

NOTE 1 Volumes of buffers can be adjusted according to protein starting amounts.

Lysis temperature should be between 60-95 °C. Visit our FAQ website for more information: www.preomics.com/faq.

SP - Storage Point: At this point, close the peptide containing tube or CARTRIDGE (use a CAP for bottom).

Peptides can be frozen at -20 °C. Storage of peptides should not exceed two weeks at -20 °C.

For extended storage, finish the protocol and store at -80 °C.

DATA ANALYSIS

Consider the following as fixed modifications in your database search:

| MODIFICATION | DESCRIPTION | COMPOSITION | SPECIFICITY | MASS | UNIMOD # |
|--------------|-----------------------------|----------------------------------|-------------|--------|----------|
| ALKYLATION | Carbamidomethyl on cysteine | C ₂ H ₃ NO | [C] | +57 Da | 4 |

Please refer to www.preomics.com for Limited Use Label License and Product Warranty.